


## **ATS 40**

### **Transfusion Set**

### **with 40 µm – Filter**




- **SAFE**
- **EASY**
- **RECOMMENDED**

## ATS 40

### Transfusion Set

ATS 40 is a high quality Transfusion product, manufactured in Germany with long standing experience in producing sterile products for hospitals. The ATS 40 has a 40µm filter and a flexible drip chamber allowing to "pump" the blood in the line.



- 40 µm Filter
- Flexible drip chamber
- 150 cm length of tubing
- Easy Roller clamp
- Luerlock-Connection

CE 0482

**Usage:** blood transfusion, specific for autologoustransfusion

List Nr.	Description	units/carton
210-000	<b>ATS-40</b> Transfusion Set with 40 µm -Filter 150 cm, Roller-clamp, Luer Lock	60