

GENESISTM RAPIDWELDTM STW

- Waferless Sterile Welding
- Touch Screen Interface
- Infrared Heating
- Unique Patented Clamps



The GENESIS RapidWeldTM STW

- yields accurate, reliable, sterile welds without the requirement of a heated consumable copper wafer
- can improve operator efficiency and maximize the shelf life of valuable blood products
- is an approved sterile connection device with the ability to produce non-contaminating welds
- is easy to use, cost effective and an eco-friendly device designed to produce sterile welds using an innovative approach to sterile tube welding



GENESIS BPS EC REP

Genesis BPS International sàrl
CP 2361 Nyon 2, Switzerland
Phone: + 41.22.363.18.16
Fax: + 41.22.363.18.17
international@genesibps.com

STRADIS-MED GmbH
Boschring 12
D-63329 Egelsbach
Amtsgericht Offenbach HRB 42891
Phone: +49 (0) 6103 73333-0
Fax: +49 (0) 6103 73333-10
E-Mail: info@stradis-med.de
Web: www.stradis-med.de

The GENESIS RapidWeld™ STW

Create strong reliable welds in 14-15 seconds • Sterile Connection Advanced Technology • Modern Design
Waferless Sterile Welding • Touch Screen Interface • Infrared Heating • Unique Patented Clamps

The RapidWeld™ STW

- yields accurate, reliable, sterile welds without the requirement of a heated consumable copper wafer
- can improve operator efficiency and maximize the shelf life of valuable blood products
- is an approved sterile connection device with the ability to produce non-contaminating welds
- is easy to use, cost effective and an eco-friendly device designed to produce sterile welds using an innovative approach to sterile tube welding

Features and Benefits

Advanced Blade Module

The stainless steel blade permits continuous and uninterrupted welding with no jamming.

- Increase operator efficiency and speed
- Notifies user when blade replacement is required
- Blade module kit: 1000 welds

User Interface with LCD Touch Screen Display

LCD color touch screen (4" x 3") display provides operator with step-by-step instructions for performing welds.

- Easy to use menu driven interface with audio & visual alarms
- Languages: English, French, German, Italian, Portuguese & Spanish

Infrared Heating Element

- Designed for continuous operation
- Maintains a high temperature to assure weld sterility

Lid Lock Mechanism

- Protects operator during welding process
- Avoids unintentional interruption of the weld

Patented Tube Holder Clamps


Unique clamp design strips fluid away from the weld site and keeps tubing securely in place.

This process can help eliminate blood cell necrosis in blood products.

- Parallel tube loading for quick and simple use
- Clearly marked indicators for operator guidance

Fast and easy step-by-step welding process

CONNECT ANY COMBINATION OF WET AND DRY PVC TUBING



Applications include, but are not limited to:

- Pathogen reduction processing - Component aliquoting - Quality control sampling - Bacterial testing - Leukoreduction - Collection set modification
- Cell washing and freezing - Apheresis set modification - Transferring research samples - Component pooling - Cellular Therapy - Other PVC applications

Article	Model	Article Code	
Genesis RapidWeld™	GRW-430	480-RW-2607-00	- Warranty: 1 year limited warranty - FDA and CE proofed
Blade Modul Kit	1000 Welds	480-RW-2478-03	

